



Standard Flame Rollator with 6-inch Wheels

The **ProBasics Flame Rollator** represents both flair & function for those who need a little more personality from their walking aid. With a stylish and affordable design fitting the majority of your patients' fundamental needs, it offers great maneuverability and stability. Constructed of lightweight aluminum and equipped with both a removable flip-back padded backrest and built-in strap for easy foldability, this rollator is designed to accompany and take your patient further whenever on the move. 300 lb weight capacity.



Water-Resistant
Nylon Pouch



Easy-to-Operate
Ergonomic Hand Brakes



6-inch Wheels

FEATURES

- Lightweight aluminum frame in burgundy and blue flame finish
- Two-button removable flip-back curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable handles for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 6" black casters with soft grip tires, ideal for indoor use
- HCPCS* Code: E0143
E0156 (seat)

COLORS

Burgundy
Flame Finish



Blue
Flame Finish



SPECIFICATIONS

Overall Width:	23.75"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	14" (W) x 12" (D)
Seat-to-Floor Height:	21"
Weight Capacity:	300 lb
Adjustable Handle Height:	30" – 36"
Wheel Size:	6"
Product Weight:	13 lb
Carton Dimensions:	24.2" x 7" x 32.6"
Shipping Weight:	16 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
RLA6BGF	ProBasics Standard Rollator with 6" Wheels	Burgundy Flame	815067072545	1/cs
RLA6BLF	ProBasics Standard Rollator with 6-inch Wheels	Blue Flame	815067072552	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.