



Junior Rollator with 6-inch Wheels

The **ProBasics Junior Rollator** is designed to fit individuals of shorter stature with a focus on style and affordability. Constructed of lightweight aluminum, it offers great maneuverability and stability. Equipped with easy-to-operate ergonomic hand brakes, a removable padded backrest and a water-resistant nylon pouch. Easily folds for storage and transport with built-in fold strap. Supports patient weights up to 250 pounds.



Water-Resistant Nylon Pouch



Easy-to-Operate Ergonomic Hand Brakes



6-inch Wheels

FEATURES

- Lightweight aluminum frame with burgundy or blue finish
- Two-button removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 6" black casters with soft grip tires, ideal for both indoor and outdoor use
- HCPCS* Code: E0143
E0156 (seat)

COLORS

Burgundy
Finish



Blue
Finish



SPECIFICATIONS

Overall Width:	22"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	13" (W) x 13.5" (D)
Seat-to-Floor Height:	19"
Weight Capacity:	250 lb
Adjustable Handle Height:	29" – 34"
Wheel Size:	6"
Product Weight:	12 lb
Carton Dimensions:	22.6" x 7.3" x 31.3"
Shipping Weight:	14.3 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
RLAJ6BG	ProBasics Junior Rollator with 6-inch Wheels	Burgundy	815067071722	1/cs
RLAJ6BL	ProBasics Junior Rollator with 6-inch Wheels	Blue	815067071715	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.